

**BIRTH DEFECTS, MYELOMENINGOCELE  
AND OMPHALOCELE  
LF01407U**



**Cleaning Instructions**

- 1.** To ensure life of product, wash with warm soapy water and completely air dry before returning to bag.
- 2.** Defect material may stain manikin; contact with ink from printed material and ballpoint pen will also stain manikin.
- 3.** Defects are to be used in simulation only.



When infants are first born, there can be many things that may go wrong. These moulage conditions will help in training a student or new parent about each of these real-life situations.

Omphalocele is when an infant's intestine or other abdominal organs are outside of the body.

Myelomeningocele (Spina Bifida) is a common birth defect. There may be neurological complications in the infant if this condition occurs.

## Instructions

To be used with the **Life/form**® Newborn Nursing Skills and ALS Simulator.

1. Place myelomeningocele into hole on lower back of infant trainer. Lubricant may be needed.
2. Place omphalocele into hole on the front of the stomach of the infant trainer. Lubricant may be needed.

