


CPR MOUTHPIECE

CLEANING & DISINFECTING INSTRUCTIONS



Instructions for cleaning and disinfecting the CPR Mouthpieces, item numbers, 2023 (all versions), 2601 (all versions), 2223, 2082, 2162, 2068, 2083, 2162B, and 2069.

Clean & Disinfect

Steps:

1. Clean as usual with your local CDC protocol.
2. To remove from solution:
 - a. hold mouthpiece with the back of it towards the palm of you hand.
 - b. squeeze the cheeks.
 - c. pour cleaning solution from the hole in the mouth.
 - d. continue squeezing the mouthpiece and shake it vigorously 2 or 3 times to assist in removing excess moisture.
3. Place mouthpiece onto its cheek on a flat surface.
4. Let air dry 1 hour to ensure that the interior of the mouthpiece thoroughly dries.



Part Numbers

2271	Brad Jr.
2275	Brad Jr. with Electronics
2801	Brad with Carry Bag
2804	Jaw Thrust Brad with Carry Bag
2805	Airway/Lung System, 24 Pack
2023	Mouth/Nose Piece, 10 Pack
2526	Carry Bag with Kneeling Pads
2280	Airway/Lung System, 24 Pack
2023ETHNIC	Mouth/Nose Piece, 10 Pack
2023EXPORT	Mouth/Lung Pieces (nares open) Export, 10 pack

Return Policy

Should it be necessary to return an item for credit or for any other reason, contact our Customer Care Department to obtain an RGA Number. Please refer to your invoice number when phoning in your request for returning merchandise. Should you have any questions or wish further information on any product we manufacture, call or write our Customer Care Department.

Warranty

Simulaids® warrants their products to be free from defects in materials and/or workmanship for a period of three years from the date of purchase, as evidenced by the date on the invoice of the product shipment to the end user. The warranty expressly does not cover abuse, accidental or purposeful damage, or any form of modification to the product. Simulaids® reserves the right to either replace affected parts or the entire unit, at their sole discretion, after investigating and reviewing the actual product and the damage. In most instances, a digital photo of the product in question showing the damage will help qualify a product for return to the factory. At no time will any product be accepted at the factory without proper return authorization issued by Simulaids®.

Freight and shipping charges are the sole responsibility of the end user. No product will be received with shipping charges due. Any product considered for warranty work must be identified by serial number and invoice number from the agency through whom the product was purchased. Without this information the product will not receive a return authorization number as required.

Do not allow print material to come in contact with the skin of your product. The ink will embed in the porous surface and be impossible to remove. Ink stains are not covered under warranty.



ONLINE:

E-MAIL:

PHONE:

ADDRESS:

NascoHealthcare.com

info@simulaids.com

1.800.431.4310

P.O. Box 1289, 16 Simulaids Dr.

Saugerties, NY 12477