

A Adam, Rouilly

AR2000 – AR2000-B

Superior Male Bedford Doll

Product Leaflet



Superior Male Bedford Doll

AR2000

AR2000-B



The Superior Male Bedford Doll is a versatile nursing manikin designed for realistic, practical clinical training. Featuring articulated joints, soft washable skin, and interchangeable anatomy, it supports procedures from patient handling and dressing to catheterisation, tracheostomy, injections, stoma care, and feeding, offering durable, hands-on practice for comprehensive nursing education.

The best-selling Superior Bedford Doll offers a range of features for both basic and advanced nursing procedures and scenarios.

From basic patient handling, dressing and care to advanced intramuscular injection, tracheostomy, catheterisation and enteral feeding the Superior Bedford Doll is an ideal full body patient manikin.

A unique soft, palpable and washable skin to ensure realism and excellent durability – ensuring years of reliable service.

Adam,Rouilly's Bedford Nurse Training Dolls offer a range of features including soft, pliable skin which can be washed, and for added realism incorporates continuous skin in elbow and knee joints.

Skills

- Articulated waist for greater flexibility
- Moulded colostomy stoma
- Bed positioning
- Patient handling
- Aspiration Bandaging – wet and dry
- Enema administration
- Intramuscular injections
- Nasogastric feeding
- Stoma care
- Stomach irrigation
- Suprapubic catheterisation
- Tracheostomy – dressing and nursing care
- Urethral catheterisation
- Washing

Immersive Learning Made Easy

Enhance the Bedford Doll with real-time vital signs, scenario control, auscultation and imaging libraries with **SimClever™**: a portable, modular “simulation-suite in a box”.

The system upgrades manikins to mid-fidelity scenario capability, all without subscription fees. Ideal for budget-conscious educators seeking immersive, flexible simulation in any clinical setting. See **Page 63** for more information.




- 1 Nasogastric Feeding
- 2 Tracheostomy Care and Dressing

Features

- Ears for aural hygiene
- Mouth and nose for nasogastric feeding, aspiration and stomach irrigation
- Articulated waist allowing body to bend and twist
- Suprapubic opening for catheterisation
- Penis for catheterisation
- Rectum for enemas
- Full range of movement through neck, shoulder, elbow, wrist, hip, knee and ankle joints
- Washable nylon rooted hair
- Moulded colostomy stoma for stoma care
- Soft feet with separate flexible toes for foot care and washing
- Replaceable intramuscular injection pad (thigh + buttock)
- Replaceable intramuscular injection pad
- Tracheostomy for draining, cleaning and dressing
- Eyes for washing and administering eye drops
- Stomach reservoir for nasogastric feeding
- Available in light and dark skin tones

Includes

Male Genitalia Only
Instruction Manual



Safety in Training, Safety in Practice

Adam,Rouilly Limited, Castle Road, Eurolink Business Park,
Sittingbourne, Kent, ME10 3AG, United Kingdom

UK

+44 (0)1795 471 378

uksales@adam-rouilly.co.uk

Export

+44 (0)1795 411 064

sales@adam-rouilly.co.uk

Adam,Rouilly has over 100 years of experience in providing quality medical models and simulators.

As part of our policy of continuous product development, the specifications of our products may change without prior notice.

Should you require further information, please contact our Sales Department, who will be pleased to assist you with your enquiry.

www.adam-rouilly.co.uk



©Copyright 2026

L2000 L1