

A Adam,Rouilly

AR305

Nose Bleed / Epistaxis Trainer

Product Leaflet



Nose Bleed / Epistaxis Trainer

AR305



The Nose Bleed/Epistaxis Trainer simulates severe nasal haemorrhages for training.

Developed with ENT specialists, it allows trainees to practise suction, nasal packing, and cautery with feedback via a light panel. Featuring variable artificial blood flow, it is an essential training tool for ENT and emergency care departments.

Skills

- Developed with the ENT Department, Addenbrooke's Hospital, Cambridge
- Interactive model to train the management of severe nosebleeds
- Highly recommended for all ENT and leading Emergency departments
- Suitable for training in cautery and nasal packing
- Cautery success at four septal sites shown by light panel (battery included)
- Trainees must locate the bleeding point, clear it by suction, then treat by packing or cautery

Features

- Capillary tubes deliver artificial blood to four septal sites at variable rates

Includes

10ml Syringes (x4)

50ml Syringe (x1)

Blood Concentrate to Make Up 2.5 Litres of Artificial Blood

Control Box and 9v Battery

Hard Plastic Head with Soft Nasal Area