

A Adam, Rouilly

AR331

The Hungry Manikin®

Paediatric Nasogastric Feeding Trainer

Product Leaflet



The Hungry Manikin® – Paediatric Nasogastric Feeding Trainer



AR331



The Hungry Manikin® Paediatric Trainer helps parents, carers, and healthcare staff safely practise nasogastric and gastrostomy feeding in children.

With transparent anatomy, fluid aspiration, pH testing, and visible tube placement, it builds confidence in measuring, inserting, and managing feeding tubes, ensuring safer care and enhanced understanding of paediatric enteral feeding.

This model was designed by nutrition nurses at Birmingham Children's Hospital to demonstrate to parents/carers the technique for inserting nasogastric and gastrostomy tubes in children.

Note

Size 8 Fr Nasogastric and Size 12 Fg Ryles tubes should be used with this model.

Skills

- Easy-to-use simulator helps identify placement stages, ensuring safe and correct passage of a nasogastric device
- Child manikin enables parents/carers to visualise tube positioning and practise related techniques
- Families report The Hungry Manikin provided invaluable hands-on practice, boosting confidence and understanding.
- Measurement, placement and securement of nasogastric tube
- Insertion of gastrostomy tube
- As the manikin accepts fluids, training in the following procedures is possible:
 - Stomach aspiration to confirm tube placement
 - Measuring pH of stomach aspirate
 - Enteral feeding via nasogastric / gastrostomy tubes
 - Administration of liquid medicines

Features

- Allows fluids to be introduced into the stomach
- Removable stomach with simple connector system for easy cleaning
- Movable epiglottis to demonstrate swallowing
- Transparent body for observing tube passage
- Divided nasal septum
- Gastrostomy opening
- Representation of lungs and oesophagus



- 1 Movement of Epiglottis
- 2 Aspiration of Stomach
- 3 Securing Nasogastric Tape

- 4 Bolus Feeding
- 5 Inflating Gastrostomy Tube
- 6 Administering via Gastrostomy

Includes

Rigid Carrying Case
Instruction Manual

Spares & Repairs

S331/8 Stomach Unit

Safety in Training, Safety in Practice

Adam,Rouilly Limited, Castle Road, Eurolink Business Park,
Sittingbourne, Kent, ME10 3AG, United Kingdom

UK

+44 (0)1795 471 378

uksales@adam-rouilly.co.uk

Export

+44 (0)1795 411 064

sales@adam-rouilly.co.uk

Adam,Rouilly has over 100 years of experience in providing quality medical models and simulators.

As part of our policy of continuous product development, the specifications of our products may change without prior notice.

Should you require further information, please contact our Sales Department, who will be pleased to assist you with your enquiry.

www.adam-rouilly.co.uk



©Copyright 2026

L0331 LI