

AirSim Advance X



[CLICK HERE TO VIEW
THE PRODUCT VIDEO](#)

WHAT IS THE AIRSIM ADVANCE X?

AirSim Advance X builds on the original AirSim design, offering enhanced capability for foundational airway training. Alongside standard procedures, it enables nasogastric (NG) tube insertion and features modular components for straightforward in-house replacement. With anatomically accurate nasal structures, an inflatable tongue, and visual feedback for lung and stomach inflation, it delivers a highly realistic training experience. AirSim Advance X includes a 3-year airway guarantee and modular design - providing long-term durability, flexibility, and excellent value for intermediate-level instruction.

WHY PURCHASE THIS MODEL?

Realistic Anatomy:

Based on CT DICOM data, with authentic anatomical landmarks (epiglottis, cricoid cartilage, tracheal rings, nasal cavity, inflatable tongue bulb).

Versatile Training:

Supports multiple procedures, including intubation, NG tube insertion, BVM ventilation and more.

Real-Time Feedback:


Visual indicators of stomach inflation and lung expansion for immediate training feedback.

Durability and Cost-Effective:

The AirSim Advance X is certified to practice 80,000+ intubation cycles without fail and is supported by a 3 year warranty.

Modular Design:

Easily replaceable parts and portable with a durable carry case.



Product Code: AA91100X
Dimensions: 47x27x34cms
Weight: 5kg

PROCEDURE TRAINING:

- Double Nasotracheal Intubation
- Bag-Valve Mask (BVM) Ventilation
- Supraglottic Device Insertion
- Direct & Video Laryngoscopy
- Awake Fibre Optic Examination
- Endotracheal Tube Insertion

IDEAL FOR TRAINING:

- Anaesthesiologists
- Emergency Medicine Residents
- Paramedics
- EMTs

UPGRADE YOUR MANIKIN:

Expand airway training with the Bronchial Tree Upgrade - an anatomically accurate fourth-generation bronchi model enabling diagnostic bronchoscopy, lung isolation, and fibreoptic intubation. Exceptional internal detail supports realistic, repeatable practice.

Product Code: BR100



TRAINING YOU CAN TRUST 3 YEAR PROMISE!

We are so confident in the quality and durability of our products that we provide a 3 year warranty on the AirSim Advance X

Product Code: AA91100X
Dimensions: 47x27x34cms
Weight: 5kg