

AirSim Bronchi



[CLICK HERE TO VIEW
THE PRODUCT VIDEO](#)

WHAT IS THE AIRSIM BRONCHI?

AirSim Bronchi is a durable, anatomically accurate airway trainer designed for early-stage and foundational skills development. Featuring a true-to-life nasal cavity and inflatable tongue, it supports realistic training in intubation, bronchoscopy, and ventilation techniques. The model provides visual feedback through lung and stomach inflation and can withstand over 80,000 intubation cycles. Ideal for beginner learners seeking a reliable and user-friendly training manikin.

WHY PURCHASE THIS MODEL?

Realistic Anatomy:

Anatomically correct nasal cavity, which allows for training in nasotracheal intubation and nasally directed fibre optic examination.

Versatile Training:

Features bronchi which provides internal anatomically correct features down to the fourth generation. Inflatable tongue bulb to create tongue edema providing a variance during training. Novel neck construction which can be secured into various positions

Real-Time Feedback:

Visual indicators of stomach inflation and lung expansion for immediate training feedback.

Durability and Cost-Effective:

The AirSim Bronchi is certified to practice 80,000+ intubation cycles without fail and is supported by a 3 year warranty.

Modular Design:

Delivered ready to use in a durable carry case for easy transportation and safe storage.

PROCEDURE TRAINING:

- Diagnostic Bronchoscopy
- Single Nasotracheal Intubation
- Bag-Valve Mask (BVM) Ventilation Techniques
- Full range of Supraglottic Device Insertion
- Direct Laryngoscopy
- Awake Fibre Optic Examination
- Endotracheal Tube Insertion
- Combi Tube Insertion
- Single Lung Isolation Techniques

IDEAL FOR TRAINING:

- Paramedics
- Anesthesiologists
- Emergency medical technicians
- Respiratory therapists

TRAINING YOU CAN TRUST 3 YEAR PROMISE!

We are so confident in the quality and durability of our products that we provide a 3 year warranty on the AirSim Bronchi.

Product Code: AA51100
Dimensions: 47 x 27 x 34 cm
Weight: 3kg