

AirSim Child Combo X



[CLICK HERE TO VIEW
THE PRODUCT VIDEO](#)

WHAT IS THE AIRSIM CHILD COMBO X?

AirSim Child Combo X is a pediatric airway and front-of-neck access trainer based on a 6-year-old anatomy. It supports intubation, ventilation, NG tube placement, and surgical airway procedures. Featuring a realistic larynx, inflatable tongue, and quick-change consumables, it's ideal for high-volume training. Certified for 80,000+ intubation cycles and backed by a 3-year warranty, it offers durability, realism, and portability.

WHY PURCHASE THIS MODEL?

Realistic Anatomy:

Based on a 6-year old. Realistic larynx which features palpable cricoid landmarks, laryngeal cartilages and tracheal rings to provide positive user feedback and easy identification of the sternal notch.

Versatile Training:

ideal for training in intubation, needle & surgical cricothyroidotomy, percutaneous tracheostomy

Real-Time Feedback:

Cost-effective consumables which are quick and easy to interchange, perfect for a busy teaching environment. Inflatable tongue bulb to create tongue edema providing a variance during training. User feedback through the visualization of stomach inflation and lung expansion.

Durability and Cost-Effective:

The AirSim Child Combo X is certified to practice 80,000+ intubation cycles without fail and is supported by a 3 year warranty.

Modular Design:

Delivered ready to use in a durable carry case for easy transportation and safe storage.

PROCEDURE TRAINING:

- Diagnostic Bronchoscopy
- Endotracheal/ Nasotracheal Intubation
- Bag-Valve Mask (BVM) Ventilation Techniques
- Full range of Supraglottic Device Insertion
- Direct Laryngoscopy
- Awake Fibre Optic Examination
- Endotracheal Tube Insertion
- Lung Isolation Techniques

Product Code: CC10006X
Dimensions: 41 x 26 x 23cm
Weight: 3kg

RECOMMENDED ADD-ON:

- Super Saver Cric Training Pack (100 Procedures) Product Code: ATP0100



ADDITIONAL INFORMATION

- Each larynx insert (NLX050) facilitates 1 cricothyroidotomy and 2-3 tracheostomy procedures.
- Each wrap around neck skin (RSN2005) facilitates 15-20 incisions.


IDEAL FOR TRAINING:

- Pediatric anesthesiologists
- Emergency medicine doctors
- Airway course instructors

UPGRADE YOUR MANIKIN:

Expand airway training with the Bronchial Tree Upgrade - an anatomically accurate fourth-generation bronchi model enabling diagnostic bronchoscopy, lung isolation, and fiberoptic intubation. Exceptional internal detail supports realistic, repeatable practice.

Product Code: BR100



TRAINING YOU CAN TRUST 3 YEAR PROMISE!

We are so confident in the quality and durability of our products that we provide a 3 year warranty on the AirSim Child Combo X

Product Code: CTC91100X
Dimensions: 41 x 26 x 23cm
Weight: 3kg