

AirSim Child X



[CLICK HERE TO VIEW
THE PRODUCT VIDEO](#)

WHAT IS THE AIRSIM CHILD X?

AirSim Child X is a realistic pediatric airway management trainer, anatomically modeled on a 6-year-old to reflect accurate proportions and dental features. Designed for high-repetition training, it supports over 80,000 intubation cycles and is backed by a 3-year airway guarantee. The model includes an inflatable tongue to simulate edema, visual feedback for lung and stomach inflation, and a modular construction for in-house part replacement. Ideal for practicing pediatric intubation, NG tube insertion, BVM ventilation, and supraglottic airway techniques in a portable, durable format.

WHY PURCHASE THIS MODEL?

Realistic Anatomy:

Based on a 6-year old with an anatomically correct airway.

Versatile Training:

Ideal for beginners in pediatric airway management techniques. Inflatable tongue bulb to create tongue edema providing a variance during training.

Real-Time Feedback:

Visual indicators of stomach inflation and lung expansion for immediate training feedback.

Durability and Cost-Effective:

The AirSim Child X is certified to practice 80,000+ intubation cycles without fail and is supported by a 3 year warranty.

Modular Design:

Modular design allowing for all parts to be changed within your organisation (if required).



Product Code: AC10006X
Dimensions: 41 x 26 x 23cm
Weight: 3kg

PROCEDURE TRAINING:

- Double Nasotracheal Intubation
- Bag-Valve Mask (BVM) Ventilation
- Supraglottic Device Insertion
- Direct & Video Laryngoscopy
- Awake Fibre Optic Examination
- Endotracheal Tube Insertion


IDEAL FOR TRAINING:

- Pediatric anesthesiologists
- Emergency medicine doctors
- Airway course instructors

UPGRADE YOUR MANIKIN:

Expand airway training with the Bronchial Tree Upgrade - an anatomically accurate fourth-generation bronchi model enabling diagnostic bronchoscopy, lung isolation, and fiberoptic intubation. Exceptional internal detail supports realistic, repeatable practice.

Product Code: BR100



TRAINING YOU CAN TRUST 3 YEAR PROMISE!

We are so confident in the quality and durability of our products that we provide a 3 year warranty on the AirSim Child X

Product Code: AC10006X
Dimensions: 41x26x23cms
Weight: 3kg