

# TruBaby X Lite



[CLICK HERE TO VIEW  
THE PRODUCT VIDEO](#)

## WHAT IS THE TRUBABY X LITE?

TruBaby X Lite is incredibly lifelike with the appearance, weight, size and movement of a 50th percentile 5-month old infant. This model is ideal for nurses and various other paediatric medical professionals practicing Directly Observed Practical Skills (DOPS), Paediatric Basic Life Support (PBLs).

## WHY PURCHASE THIS MODEL?

### **Realistic Anatomy:**

Appearance, weight, size, and movement of a 50th percentile 5-month old infant.

### **Versatile Training:**

Train a number of nursing procedures such as the new and improve IV cannulation in the hand, arm & foot or urethral catheterisation.

### **Durability and Cost-Effective:**


All parts are cost effective and can be replaced quickly when required. Initial set-up time is less than 5 mins. The TruBaby X Lite is certified to practice 80,000+ intubation cycles without fail and is supported by a 2 year warranty.

### **Portable:**

Delivered in a durable carry case for easy transportation and safe storage.

## PACKAGE CONTENTS

- x4 sets of IO inserts
- x1 male genitalia insert (female insert attached when delivered)
- x1 bottle of TruCorp lubrication (100ml)
- x1 bottle of artificial blood concentrate (250ml)



Product Code: TBL10001X  
Dimensions: 79 x 46 x 28cm  
Weight: 8.5kg

## IDEAL FOR TRAINING:

- Airway Management techniques (oral and nasal)
- Peripheral Venous Cannulation (hand, arm and foot)
- IO Tibia
- Intramuscular Injection
- Urethral Catheterization (male and female)
- Cardiopulmonary Resuscitation (CPR)

## FEATURES & BENEFITS:

- IO Tibia: Realistic resistance when penetrating the medullary cavity. Single use insert provides each learner with a lifelike training experience
- Urethral Catheterization: Interchangeable male and female inserts to practise both gender techniques. Upon successful entry, fluid will flow from the catheter
- CPR: Full recommended depth of 1.5 inches can be achieved with full chest re-coil



### TRAINING YOU CAN TRUST- 2 YEAR PROMISE!

We are so confident in the quality and durability of our products that we provide a 2 year warranty on the TruBaby X Lite.

Product Code: TBL10001X  
Dimensions: 79 x 46 x 28cm  
Weight: 8.5kg