

TruInfant IO Leg & Consumable Bundle



[CLICK HERE TO VIEW THE PRODUCT VIDEO](#)

WHAT IS THE TRUINFANT IO LEG?

The TruInfant IO Leg is a highly realistic task trainer designed for healthcare professionals to practice and perfect critical pediatric skills. This model is based on a 5-month-old infant 50th percentile, it provides an authentic training experience in three main areas: Intraosseous (IO) infusion (tibia), IV cannulation (foot), and intramuscular injections (thigh). The TruInfant IO Leg comes with a fully assembled IO Leg and a consumable bundle that includes 40 IO Inserts.

WHY PURCHASE THIS MODEL?

Versatile Training:

TruInfant IO Leg is perfect for training as it allows for 3 procedures in one unit. Replaceable legs to extend the trainer's life and maintain training quality.

Real-Time Feedback:

Realistic blood flashback when the procedure is performed successfully.

Durability and Cost-Effective:

Quick and easy replacement of IO inserts. Each IO insert can withstand up to 5 procedures. 250+ needle penetrations per IV and IM injection site.

Quick & Easy Set-Up:

Ready for use in under 5 minutes.

PROCEDURE TRAINING:

- IO tibia needle placement.
- IV cannulation in the foot.
- Intra-muscular injection sit

IDEAL FOR TRAINING:

- Paramedics
- Pediatricians
- Pediatric Nurses
- Emergency medicine professionals

TRAINING YOU CAN TRUST 1 YEAR PROMISE!

We are so confident in the quality and durability of our products that we provide a 1 year warranty on the TruInfant IO Leg

Product Code: TINIOLEG2
Dimensions: 41 x 26 x 23 cm
Weight: 7.5kg

Skill-Specific Training:

Intraosseous (IO) Tibia Needle Placement:

- **Realistic Anatomy:** Clearly defined tibial tuberosity and patella guide proper needle placement.
- **Lifelike Resistance:** Penetrating the medullary cavity mimics the feel of real bone tissue.
- **Positive Feedback:** Simulated blood flow confirms correct needle placement and technique.



IV Cannulation (Foot):

- **Comprehensive Venous System:** Includes dorsal venous arch, great saphenous vein, and lesser saphenous vein.
- **Realistic Feel:** Provides lifelike resistance during cannulation and blood flashback.
- **Multiple Cannulation Sites:** Allows practice on various venous locations for comprehensive training. Allowing up to 250+ needle penetrations.



Intramuscular Injections (Thigh):

- **Anatomical Accuracy:** Ensures correct identification and technique for intramuscular injections.
- **Realistic Resistance:** Simulates the feel of injecting into muscle tissue and can withstand up to 250+ needle penetrations.



Product Code: TINIOLEG2
Dimensions: 41 x 26 x 23 cm
Weight: 7.5kg