

TruTourniquet



[CLICK HERE TO VIEW
THE PRODUCT VIDEO](#)

WHAT IS THE TRUTOURNIQUET?

The TruTourniquet has been designed to practise massive haemorrhage control in lower extremity trauma. The model is based on a right thigh and provides realistic bleeding control with visual feedback on user performance.

WHY PURCHASE THIS MODEL?

Realistic Anatomy:

TruTourniquet is a realistic and robust above-knee blast injury task trainer. It features one bleeding point and a fractured femur.

Versatile Training:

Users can simulate a massive haemorrhage when the product is connected to the pressurized blood pump. Suitable for single or double tourniquet application.

Real-Time Feedback:

The model provides visual feedback when the correct tourniquet force has been applied.

Durability and Cost-Effective:

Quick & easy set-up - less than 5 minutes. If the blood system is empty simply pour the blood back into the pump.

Suitable for All Training Environments:

Realistic weight and vertical stump movement to practise full application procedure. Ability to detach stump from the base to use in a non-classroom environment. The splash shield helps avoid blood overflow and keeps mess to a minimum. It can be easily attached/ detached as required.

IDEAL FOR TRAINING:

- Military Medics
- Paramedics and EMTs
- Trauma and Emergency Care Physicians
- Nurses specialising in emergency and trauma care

Tested for over 40,000+
tourniquet applications
without fail

TRAINING YOU CAN TRUST 1 YEAR PROMISE!

We are so confident in the quality and durability of our products that we provide a 1 year warranty on the TruWound.

Product Code: TTQ1000
Dimensions: 47 x 27 x 34cm
Weight: 7.5kg