

TruWound



[CLICK HERE TO VIEW
THE PRODUCT VIDEO](#)

WHAT IS THE TRUWOUND?

TruWound is a realistic gunshot wound model with a pressure blood supply system to simulate severe haemorrhaging. It allows realistic wound packing training techniques with visual feedback when the correct pressure is applied to the packed wound. The pressurized blood system allows the instructor to control blood flow to either 1, 2 or all 3 wounds at once, making it ideal for a fast paced classroom environment.

WHY PURCHASE THIS MODEL?

Realistic Anatomy:

Size, appearance, and feel is realistic to a large bullet wound in the upper thigh and can be packed with haemostatic pads to teach wound packing procedures.

Versatile Training:

Users can simulate severe bleeding with venous and arterial flows with the pressurised blood system.

Real-Time Feedback:

The model provides visual feedback when the correct pressure has been applied to the packed wound.

Durability and Cost-Effective:

Easy to set up. Initial set-up is less than 5 minutes. The TruWound blood collection containers help keep mess to a minimum during packing procedures and the artificial blood is washable and will not stain the model.

Increased Efficiency:

Suitable for a classroom setting by having three wounds operating simultaneously.

IDEAL FOR TRAINING:

- Military Medics
- Paramedics and EMTs
- Trauma and Emergency Care Physicians
- Nurses specialising in emergency and trauma care

It's durable!
Tested for over 1,500+
procedures without fail

TRAINING YOU CAN TRUST 1 YEAR PROMISE!

We are so confident in the quality and durability of our products that we provide a 1 year warranty on the TruWound.

Product Code: TWO1000
Dimensions: 47 x 27 x 34cm
Weight: 4.5kg